



Making it Easier
to Handle
**Compression
Therapy**



Application Aids and Accessories
for JOBST® Compression Garments

Arion Application Aids

JOBST® is the preferred distributor of Arion products.

Arion products are made from very smooth material, treated with a special coating, which together with double-layer product design create an extremely smooth gliding product, thus significantly minimizing the friction between the compression garment and skin.

For putting ON and taking OFF compression garments


Arion Magnide® 2in1 On & Off

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight



Size	Shoe size	BNR	UOM / Box
Medium	< 7,5	7957110	1
Large	≥ 8 - 10	7957111	1
X-Large*	≥ 10,5	7957112	1

* XL also suitable for extra wide feet with cA up to 28 cm




Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product



Size	Shoe size	BNR	UOM / Box
Small	< 3 - 6	7063301	1
Medium	6.5 - 8	7063302	1
Large	8.5 - 10.5	7063303	1
X-Large	≥ 11	7063304	1



Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®



Size	Circumference of lower leg	BNR	UOM / Box
Medium	≤ 40 cm (15.7")	7957106	1
Large	≤ 60 cm (23.6")	7957107	1
Extra large	≤ 100 cm (39.5")	7957107	1





For putting ON compression garments

Arion Easy-Slide on

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application



Size	Shoe size	BNR	UOM / Box
X-Small	≤ 2	7965803	1
Small	2.5 - 5.5	7965804	1
Medium	6 - 8	7965802	1
Large	8.5 - 11	7965902	1
X-Large	≥ 11.5	7966001	1





Arion Easy-Slide Arm on

- The user-friendly application aid makes putting on compression arm sleeves quick and easy
- A straightforward donning method in combination with the application aid



Size	Circumference of widest part of the arm	BNR	UOM / Box
Medium	37 - 38.5 cm (14.5" - 15.1")	7966102	1
Large	39 - 41 cm (15.3" - 16.1")	7510001	1



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes



Size	Shoe size	BNR	UOM / Box
Medium	≤ 7.5	7957100	1
Large	8 - 12	7957101	1
X-Large*	≥ 13	7957102	1




Arion Easy-Slide

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children



Size	Age	Length	BNR	UOM / Box
X-Small	0-3 yrs	29.5 cm (11.6")	7957108	1
Small	3-6 yrs	39.5 cm (15.5")	7957109	1




For putting on above-knee prosthesis with suction socket

Arion Proth-Pro


- Specifically designed to help alleviate the problems associated with putting on a prosthesis
- One single seam results in a very strong product with an improved lifetime
- Light, smooth material minimizes friction and reduces the pain caused by friction when putting on the prosthesis



Size	Circumference		BNR	UOM / Box
	Proximal	Distal		
X-Small	< 41 cm (16.1")	25 cm (9.8")	6032800	1
Small	41-45 cm (16.1"-17.7")	28 cm (11")	6032801	1
Medium	46-58 cm (18.1"-22.8")	30 cm (11.8")	6032802	1
Large	59-70 cm (23.2"-27.6")	32 cm (12.6")	6032803	1
X-Large	71-80 cm (28"-31.5")	34 cm (13.4")	6032804	1
XX-Large	> 81 cm (31.5")	36 cm (14.2")	6032805	1
Extender	One size fits all - 120 cm (47.2")		6032806	1



How to use the Extender



JOBST® Donning Devices and Accessories

Help your patients put on JOBST® garments, while providing convenience and value.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.



Size	BNR	UOM / Box
Small	131202	1 pair
Medium	131203	1 pair



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles

Size	BNR	UOM / Box
Standard	110913	1 donner



Warnings for Arion Application Aids

- Do not step on the product. The highly slippery material makes it a slipping hazard.
- The Magnide[®], Magnide[®] 2in1 and the SlideX[®] contain several magnets. For this reason, keep this product away from pacemakers and defibrillators. The magnetic field may interfere with the proper operation of this equipment (minimum distance = 30 cm)
- When using the product it is important to use the correct size for you. Using an incorrect size will prevent the product from working properly, and it may become damaged or wear out too quickly. The correct size is determined by your shoe size.
- Wear special rubber gloves when donning or doffing. This gives you maximum grip on the stocking, and reduces the physical strength needed.
- This product is guaranteed for a period of six months from the date of purchase against defects resulting from faulty material or construction. Within the guarantee period, the product will be exchanged free of charge in the case of manufacturing and/or material defects. Any exchange under the guarantee will only be carried out if it can be demonstrated by means of a purchase receipt that the date on which the complaint was filed falls within the guarantee period. To use this service, please contact your supplier. You are asked to return the product washed and to follow the washing instructions.

 @jobst.canada

 jobstcanada.com



JOBST[®],
an Essity brand

136167 63888E R1 ©Essity Hygiene and Health AB, 2020 H20

Essity
1275 North Service Road West, Suite 800
Oakville, ON Canada L6M 3G4
Tel. 1-877-978-5526 Fax 1-877-978-9703