



Making it Easier  
to Handle  
**Compression  
Therapy**



Application Aids and Accessories  
for JOBST® Compression Garments

# Arion Application Aids

JOBST® is the preferred distributor of Arion products.

Arion products are made from very smooth material, treated with a special coating, which together with double-layer product design create an extremely smooth gliding product, thus significantly minimizing the friction between the compression garment and skin.

## For putting ON and taking OFF compression garments


### Arion Magnide® 2in1

- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight



Size	Men's Shoe size	Women's Shoe size	BNR	UOM / Box
Medium	≤ 6	≤ 7.5	7957110	1
Large	6.5 - 10	8.5 - 11	7957111	1
X-Large*	≥ 11	-	7957112	1

\*XL also suitable for extra wide feet with a circumference up to 11" / 28 cm




### Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product



Size	Men's Shoe size	Women's Shoe size	BNR	UOM / Box
Small	3 - 4.5	≤ 4 - 6	7063301	1
Medium	5 - 6.5	6.5 - 8	7063302	1
Large	7 - 8	8.5 - 9.5	7063303	1
X-Large*	≥ 8.5	≥ 10	7063304	1





# For putting ON compression garments

## Arion Easy-Slide on

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application

Size	Men's Shoe size	Women's Shoe size	BNR	UOM / Box
X-Small	≤ 2.5	≤ 3.5	7965803	1
Small	3 - 4.5	4 - 6	7965804	1
Medium	5 - 6.5	6.5 - 8	7965802	1
Large	7 - 8	8.5 - 9.5	7965902	1
X-Large*	≥ 8.5	≥ 10	7966001	1





## Arion Easy-Slide Arm on

- The user-friendly application aid makes putting on compression arm sleeves quick and easy
- A straightforward donning method in combination with the application aid

Size	Arm Length*	Max. Circumference**	BNR	UOM / Box
Medium	23" - 26.5" (59-67cm)	18" (45.6cm)	7966102	1
Large	>26.5" (>67cm)	22.5" (57.6cm)	7510001	1

\*Arm length from the tip of the middle finger to the armpit  
 \*\*For a full-length distribution of the Easy-Slide Arm over your arm, we recommend that the circumference of the widest part of the arm is min. 0.4"/1cm smaller than the max. product circumference



## Arion Magnide® on

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes

Size	Men's Shoe size	Women's Shoe size	BNR	UOM / Box
Medium	≤ 6	≤ 7.5	7957100	1
Large	6.5 - 10	8.5 - 11	7957101	1
X-Large*	≥ 11	-	7957102	1




\* XL also suitable for extra wide feet with a circumference up to 11" / 28 cm

## Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide® (2in1)



Size	Circumference of lower leg	BNR	UOM / Box
Medium	≤ 15.7" (40cm)	7957106	1
Large	≤ 23.6" (60cm)	7957107	1
Extra large	≤ 39.4" (100cm)	7957118	1

## Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children



Size	Age	Length	BNR	UOM / Box
X-Small	0-3 yrs	11.5" (29.5cm)	7957108	1
Small	3-6 yrs	15.5" (39.5cm)	7957109	1



## For putting on above-knee prosthesis with suction socket


### Arion Proth-Pro

- Specifically designed to help alleviate the problems associated with putting on a prosthesis
- Light, smooth material minimizes friction when putting on the prosthesis



Size	Circumference		BNR	UOM / Box
	Proximal	Distal		
X-Small	15.5" (40cm)	4.5" (12cm)	6032800	1
Small	16"-17.5" (41 - 45cm)	5" (13cm)	6032801	1
Medium	18"-23" (46 - 58cm)	5" (13cm)	6032802	1
Large	23.5"-27.5" (59 - 70cm)	6" (15.5cm)	6032803	1
X-Large	28"-31.5" (71 - 80cm)	6" (16cm)	6032804	1
XX-Large	32" (81cm)	-	6032805	1
Extender	One size fits all - 47.2" (120cm)		6032806	1

### How to use the Extender



# JOBST® Donning Devices and Accessories

Help your patients put on JOBST® garments, while providing convenience and value.

## JOBST® Easy Wash & Wear Kit

- Everything needed to help apply and extend the life of JOBST® Hosiery
- Kit includes:
  - JOBST® Donning Gloves
  - Wash Bag
  - Jolastic® Washing Solution
  - Helpful hints for putting on JOBST® Hosiery brochure



Size	BNR	UOM / Box
Includes gloves, Jolastic, wash bag and brochure	112004	1

Your Guide to Healthy Legs Brochure  
PN: 61834

## JOBST® Donning Gloves\*

- Help to create friction and provide grip
- Help to evenly distribute fabric on the leg

\* Caution: This product contains natural rubber latex which may cause allergic reactions.





Size	BNR	UOM / Box
X-Small	7740027	1 pair
Small	7740028	1 pair
Medium	7740029	1 pair
Large	7740030	1 pair
X-Large	7740031	1 pair

## JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box
Standard	110913	1 donner
Large	110910	1 donner




## JOBST® Stocking Slips

- For use with open toe styles only
- Helps the stocking slide smoothly onto the foot
- Reduces friction between the patient's skin and the compression stocking



Size	BNR	UOM / Box
1-size only	7846301	1



## JOBST® Roll-On Adhesive

- Designed to keep JOBST® Hosiery in place
- Ideal for thigh high lengths that do not have a silicone band or when a garter belt is not an option
- Helpful Hint: For flat-knit garments, apply in vertical strips

Size	BNR	UOM / Box
2 fl. oz. Bottle	7927410	1
Case of 12 Bottles	7927411	12



## Jolastic® Washing Solution

- Warm water wash formulated to remove soil, body oils, and skin salts
- Acts quickly and thoroughly without damage to fabric
- Helps extend the life of JOBST® garments

Size	BNR	UOM / Box
4 fl. oz. Bottle	130999	10
12 fl. oz. Bottle	112000	1
1 qt. Bottle	112003	1



## Dycem® Non-Slip Material

- Improves grip
- The increased friction of the anti-slip material provides hold and resistance when donning a garment
- Non-slip surface for easier donning

Size	BNR	UOM / Box
One size	7763600	1



# For Arion Application Aids

## CONTRAINDICATIONS

- Do not use the product if you are allergic to one of its component.
- Open wounds should be covered with an appropriate dressing before using the product.

## WARNINGS

- The Magnide®, Magnide® 2in1 and the SlideX contains several magnets. Keep this product away from pacemakers and defibrillators. The magnetic field may interfere with the proper operation of this equipment (minimum distance = 11.8"/30 cm).
- Suffocation hazard: Keep the product away from babies, children, and pets.
- Slipping hazard: Do not step on the product. Highly slippery material.
- For external use only / Intended for skin contact. Seek medical attention if irritation develops.

Scan here to  
view the Arion  
instructional videos:



or go to: [uqr.to/JobstUS-videos](https://uqr.to/JobstUS-videos)

 /JOBSTUSA

 @JOBSTforUSA

 @JOBST\_USA

 [jobst-usa.com](https://jobst-usa.com)



JOBST®,  
an Essity brand

BSN Medical Inc., an Essity company  
5825 Carnegie Blvd., Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325  
To order toll-free: Tel. (+1) 800 537 1063

63888 R2 ©2023 BSN Medical Inc. A23