


Patient Name / ID Code or File # \_\_\_\_\_ DOB \_\_\_\_\_ Date \_\_\_\_\_

Address \_\_\_\_\_ Gender M  F

City/State/Zip \_\_\_\_\_

Document the measurements for the **stump limb only**



Left Right

### Circumferences

Use Elvarex landmarks

cG<sup>++/+\*\*\*</sup> \_\_\_\_\_  
 cF<sup>++</sup> \_\_\_\_\_  
 cE<sup>+</sup> \_\_\_\_\_  
 cD<sup>+/0\*\*\*</sup> \_\_\_\_\_  
 cC<sup>++</sup> \_\_\_\_\_  
 cB1<sup>++</sup> \_\_\_\_\_  
 cB<sup>+</sup> \_\_\_\_\_  
 cY<sup>0</sup> \_\_\_\_\_  
 cA<sup>+</sup> \_\_\_\_\_

### Lengths

Take lengths from end of stump. End of stump = 0

/G \_\_\_\_\_  
 /F \_\_\_\_\_  
 /E \_\_\_\_\_  
 /D \_\_\_\_\_  
 /C \_\_\_\_\_  
 /B1 \_\_\_\_\_  
 /B \_\_\_\_\_  
 /A (medial) \_\_\_\_\_  
 /A (lateral) \_\_\_\_\_  
 /Z \_\_\_\_\_

0 no tension  
 + light tension  
 ++ heavy tension



With plate



Without plate

### Circumference around the base of stump

This will help give the end of the garment better shape.

\_\_\_\_\_ cm

### Hemispheric Lengths of end of stump

Taken from the base of stump circumference mark

Front to Back \_\_\_\_\_ cm

Left to Right \_\_\_\_\_ cm

The illustration above is an example of where a stump garment might end. Measure all landmark circumferences and lengths until the end of your patient's limb. Document the measurements for the **stump limb only**. Measurements for the full limb should be documented on Elvarex Lower Extremity order form (50333). If a stump garment is needed for both legs, please complete two stump order forms (this form).

#### ADDITIONAL INFORMATION:

- Include a tracing of the end of the stump to show any abnormal protrusions.  
- Tracing must contain 1" x 1" scale for reference.
- Pictures are very helpful: front, side, back.
- Please have fitter call Elvarex Customer Service Team before measuring if assistance is needed:  
1-800-221-7573, option 3.

**CAUTION:** This product contains natural rubber latex which may cause allergic reactions.

- \* Design Pressure
- \*\* If measuring is done in lying position, cA please apply 0 tension
- \*\*\* If silicone band & straight ending